

Preventing Wine Spoilage:

Rapid Screening & Quantitative Analysis of Off-flavor Phenolic Compounds by DART Mass Spectrometry

Elizabeth Crawford^{1,2}, Paola Domizio^{3,4},

Brian Musselman¹, C. M. Lucy Joseph³, Linda F. Bisson³,

Bart C. Weimer⁵ and Richard Jeannotte^{5,6}

¹IonSense, Inc, Saugus, MA, USA

²Dept. of Food Analysis & Nutrition, Institute of Chemical Tech., Prague, Czech Republic

³Dept. of Viticulture & Enology, Univ. of California-Davis, Davis, CA, USA

⁴Dipart. di Gestione Sistemi Agrari, Alimentari e Forestali (GESAAF), Univ. degli Studi di Firenze, Italy

⁵Dept. of Health & Reproduction, School of Veterinary Med, Univ. of California-Davis

⁶Facultad de Ciencias, Univ. de Tarapacá, Arica, Chile

11th June 2013


2013 ASMS Annual Meeting

Minneapolis, MN, USA

Brettanomyces

Background:

- Budding yeast found widely distributed in nature
- Discovered in beer in 1904 (Claussen), in wine (Krumbholz & Tauschanoff, 1930) and again in 1940 (Custers)
- Produces a wide array of aromatic compounds
- Wine cellar contamination was widespread
- **“Brett” characters can compete with varietal characters for dominance of wine profile**



When Is It Spoilage?

- High concentration, dominating wine profile
- Conflict with wine matrix characters
- Suppression of varietal character
- Enhancement of off-notes
- Lactic acid bacteria often found in wines with *Brettanomyces*

Recovery Thresholds:

- Chatonnet* has defined spoilage as:
 - >426 ppb of 4-ethylphenol (4-EP) and 4-ethylguaiacol (4-EG)
 - >620 ppb of 4-EP
- 50% of tasters can detect 605 ppb in wine or 440 ppb in water of 4-EP

* Chatonnet, P.; Boidron, J. N.; Dubourdieu, D. Influence des conditions d' élevage et de sulfitage des vins rouges en barriques sur leur teneur en acide acétique et en éthyl-phenols. *J. Int. Sci. Vigne Vin.* **2003**, *27*, 277-298.

Main Aromas and Incidence of Spoilage


Associated Aromas:

- Band-Aid (4-Ethyl Phenol)
- Earthy (Geosmin)
- Horsy
- Leather
- Putrid
- Soy
- Tobacco

Country	>426ppb	>620ppb
France	36%	28%
Italy	49%	19%
Australia	59%	46%
Portugal	42%	27%

Wines may contain up to 50 ppm (!) of 4-EP

Evolving Open Air Ionization



Venter, A.; Nefliu, M. Cooks, R.G., "Ambient Desorption Ionization Mass Spectrometry", Trends in Anal. Chem., 27, 284-290, 2008.

Classic vs. Ambient Ionization

GC-MS Method:

Sample Preparation:

- **Total runtime = ~50 minutes**
 - **Total sample prep time = 15 minutes** with **100% analyst involvement**
 - **GC-MS runtime of 32 minutes**
 - Internal standards (4-EP-d4; 4-EG-d5) spiked into wines (5 mL), vortexed, solvent added (t-butyl ether), centrifugation, 5 μ L injection onto GC
- (Adapted from Rayne & Eggers (2008) (Am. J. Enol. Vitic. 59:92-97))

DART Ambient Ionization Method:

- **Direct sample analysis** under ambient conditions
 - **No separation** of sample components (chromatographic)
 - **Minimal** sample preparation requirements (sample concentration)
 - Analysis time of **3 minutes** per sample

AMBIENT MASS SPECTROMETRY

CLASSIC APPROACH (GC-MS)




SAMPLE




SAMPLE PREPARATION



IONIZATION & DETECTION



CHROMATOGRAM



MASS SPECTRUM


AMBIENT MS APPROACH (DART-MS)



SAMPLE



IONIZATION & DETECTION



CHROMATOGRAM



MASS SPECTRUM

Key DART Source Parameters

- **Source Parameters:**
 - Ionization gasses
 - N₂ (Standby)
 - He (Run Mode)
 - **Temperature** of ionization gas
 - Sample introduction **speed**
- **Ions typically observed in DART-MS mass spectra:**
 - [M+H]⁺, [M-H]⁻
 - M⁺, M⁻ (ionic compounds)
 - Use of dopants to promote ionization:
 - [M+NH₄]⁺ (ammonia)
 - [M+Cl]⁻ (e.g. CH₂Cl₂)
 - [M+CF₃COO]⁻ (trifluoroacetic acid)




Figure courtesy of Prof. Jana Hajšlová
ICT Prague, Czech Republic

DART-SVP Coupled with Q Exactive

Targeted MS² Settings

- **Scan Parameters:**
 - **Negative Ion Mode**
 - **Resolution:** 35,000 FWHM at m/z 200
 - **NCE Fragmentation:** 45 (4-EP); 25 (4-EG)
 - **AGC Target:** $2e^5$ charges
 - **Max IT:** 250 ms
 - **Isolation Width (m/z):** 3.0
- **Q Exactive Inlet Parameters:**
 - **Capillary Temp:** 200° C
 - **S Lens:** 50
- **External Mass Calibration**

- **All of the following parameters were set to zero:**
 - Sheath Gas Flow, Aux Gas Flow, Sweep Gas Flow, Spray Voltage

DART-SVP Settings

- **DART Source:**
 - **Negative Ion Mode**
 - **Heater Temperature Gradient:** 100 - 300° C



* SVP = Standardized Voltage and Pressure

Classic DART Sampling: Wine on Glass Tips




Spiked Wine Sampled Direct on Glass

DART_tMS2_ChiantiWine_4-EP&4-EG_0.5-1...
3 uL spot DIP-it Method; He; UHV 3.45e-9


1/25/2013 3:53:20 PM

4-EP 45NCE;4-EG 35NCE Blue 3; 300 C

RT: 0.00 - 4.26




Negative Ion Mode



Targeted MS/MS experiment looking at major fragment ions, with direct sampling from merlot wine spiked with 4-EP and 4-EG

Gerstel Twister Stir Bar: Sample Concentration



Twister was spun overnight in orange juice with **50 ppb** spiking level of **10 pesticides**

Real Wine Sample - Direct Twister



Preliminary Results - No Glass Chamber

DART_tMS2_Wine_4-EP&4-EG_Wine without...
Gerstel Twisters 24hr spin; He; UHV 3.45e-9

1/25/2013 4:20:00 PM

4-EP 45NCE;4-EG 35NCE Blue 3; 300 C

RT: 0.59 - 2.70



Signal from real wine sample could not be detected directly, need to concentrate analytes.

Wine Sampled with PDMS Twister




**Twister stir bar
(PDMS, 0.5 mm thickness)
introduced into 2 mL wine**

**Minimum 30 minutes
unattended stirring**

MAJOR IMPROVEMENT

Twister placed in glass chamber = COMPLETE desorption




Typical DART-MS Chronogram from Stir Bar

Final Results - Glass Chamber Yields Much Improved Sampling Reproducibility

25 ppb Spiking Level of Phenolics


RT: 0.00 - 5.60




Major improvement using glass chamber for even desorption from a single Twister stir bar

Twister Stir Bar: Quantitative Results

Merlot Spiked with 4-Ethylphenol (4-EP)



Merlot Spiked with 4-Ethylguaiaicol (4-EG)



DART MS/MS Method: Figures of Merit

4-EP Conc. with IS (µg/L)	1	10	25	100	250	750
1	ND	5.9	32.8	101.1	259.0	737.7
2	ND	9.8	24.4	88.4	253.4	764.9
Mean	N/A	7.8	28.6	94.8	256.2	751.3
%Bias	N/A	-21.7	14.5	-5.2	2.5	0.2
n	2	2	2	2	2	2



4-EP Conc. (µg/L)	Low	High
	50	500
1	47.7	492.2
2	58.7	509.2
3	52.1	N/A
4	41.0	N/A
Mean	49.9	500.7
%CV	14.9	2.4
%Bias	-0.3	0.1
n	4	2

Compare: Calculated Levels of 4-EP & 4-EG

Wine Sample	DART HRAM MS/MS		GC MS	
	4-EP (µg/L)	4-EG (µg/L)	4-EP (µg/L)	4-EG (µg/L)
Sample 04	854 *	197	845	203
Sample 05	518	157	563	161
Sample 06	52	ND	129	14
Sample 09	ND	ND	110	13
Sample 14	2774 *	492	2534	433

Yellow = Brettanomyces

* Levels above selected calibration range

Twister: Effect of Time on Concentration

Wine Sample	4-EP ($\mu\text{g/L}$)	4-EG ($\mu\text{g/L}$)
Sample 44 (30 min)	3775*	287
Sample 44 (1 hour)	3276*	295
Sample 44 (12 hour)	3587*	309
RSD %	7.0	3.6
Sample 19 (30 min)	2230*	450
Sample 19 (12 hour)	2342*	440
RSD %	3.5	1.6
Sample 36 (30 min)	1443*	289
Sample 36 (12 hour)	1277*	296
RSD %	8.6	1.7

* Levels above selected calibration range

Summary

- **Remarkable correlation** between classic GC-MS method and novel DART-MS/MS method
- **Effective and labor-free and solvent-free** analyte concentration onto the **Gerstel Twister sorptive stir bars**
 - Direct screening of volatile **phenolic compounds** at **low levels (10 - 50 ppb range)** with **3 minute** sample analysis time with DART ionization
- The stir bar **DART HRAM MS/MS** method is **quantitative** over the targeted concentration range with good precision and accuracy
 - **Limits of quantification** between **10 - 50 ppb**
- **Minimized analyst interaction** with the samples with DART method:
 - Reduces analyst error
 - Increases productivity and sample throughput
 - **Real-time monitoring = quick response to changes during wine fermentation to prevent wine spoilage (minimize revenue loss)**
- **Future work:** Apply this method for broader phenolic characterization of wines for rapid wine characteristic fingerprinting

Questions?

Many thanks to...

Thermo Fisher Scientific

- Dr. Catharina Crone
- Dr. Markus Kellmann
- Dr. Tabiwang Arrey



Email: crawford@ionsense.com

Recent Publications

- **Direct analysis in real time mass spectrometry and multivariate data analysis: A novel approach to rapid identification of analytical markers for quality control of traditional Chinese medicine preparation.** Zeng S, Wang L, Chen T, Wang Y, Mo H and Qu H. *Analytica Chimica Acta* **733**:38–47, 2012
- **Rapid detection of alkaloids in Ipecac by direct analysis in real time tandem mass spectrometry (DART-MS/MS).** Sun L, Hu X, Liu L, Jin H and Lin R. *Zhongguo Zhong Yao Za Zhi* **37**:1426–30, 2012 (Article in Chinese)
- **Metabolomic fingerprinting employing DART-TOF MS for authentication of tomatoes and peppers from organic and conventional farming.** Novotna H, Kmiecik O, Gałazka M, Krtkova V, Hurajová A, Schulzová V, Hallmann E, Rembiałkowska E and Hajšlova J. *Food Additives & Contaminants: Part A* **29**(9):1335–1346, 2012
- **Rapid identification of synthetic cannabinoids in herbal samples via direct analysis in real time mass spectrometry.** Musah R, Domin M, Walling M and Shepard J. *Rapid Commun. Mass Spectrom.* **26**(9):1109–1114, 2012
- **Analysis of isoflavones in soybeans employing direct analysis in real-time ionization–high-resolution mass spectrometry.** Lojza J, Cajka T, Schulzova V, Riddellova K and Hajslova J. *J. Sep. Science.* **35**(3):476–481, 2012
- **Confined direct analysis in real time ion source and its applications in analysis of volatile organic compounds of Citrus limon (lemon) and Allium cepa (onion).** Li Y. *Rapid Commun. in Mass Spectrom.* **26**(10):1194–1202, 2012
- **Chemometric Classification of Morphologically Similar Umbelliferae Medicinal Herbs by DART-TOF-MS Fingerprint.** Lee SM, Kim HJ and Jang YP. *Phytochem. Anal.* **23**(5):508-512, 2012.

<http://www.ionsense.com/pdfs/DARTFoodBook18Sep12sm.pdf>